

# Biomedical and Signal Processing Lab

## Experiment 9

### Simulation of ECG pulse missing detector

**Aim:** To Simulate ECG pulse missing detector

**Objective:**


- i. To understand various wave complexes involved in ECG analysis.
- ii. To understand the scheme for measurement of heart rate from ECG.
- iii. To understand the significance of missing QRS wave and its detection.

**Theory:**

The electrocardiogram (ECG) is a diagnostic tool that measures and records the electrical activity of the heart in detail. Being able to interpretate these details allows diagnosis of a wide range of heart problems.

### ECG Interpretation

The graph paper that the ECG records on is standardized to run at 25mm/second, and is marked at 1 second intervals on the top and bottom. The horizontal axis correlates the length of each electrical event with its duration in time. Each small block (defined by lighter lines) on the horizontal axis represents 0.04 seconds. Five small blocks (shown by heavy lines) is a large block, and represents 0.20 seconds.




Duration of a waveform, segment, or interval is determined by counting the blocks from the beginning to the end of the wave, segment, or interval.

P-Wave: represents atrial depolarization - the time necessary for an electrical impulse from the sinoatrial (SA) node to spread throughout the atrial musculature.

- Location: Precedes QRS complex  
Amplitude: Should not exceed 2 to 2.5 mm in height  
Duration: 0.06 to 0.11 seconds

P-R Interval: represents the time it takes an impulse to travel from the atria through the AV node, bundle of His, and bundle branches to the Purkinje fibres.

- Location: Extends from the beginning of the P wave to the beginning of the QRS complex  
Duration: 0.12 to 0.20 seconds.



QRS Complex: represents ventricular depolarization. The QRS complex consists of 3 waves: the Q wave, the R wave, and the S wave.


- The Q wave is always located at the beginning of the QRS complex.  
It may or may not always be present.

The R wave is always the first positive deflection.

The S wave, the negative deflection, follows the R wave

- Location: Follows the P-R interval  
Amplitude: Normal values vary with age and sex

Duration: No longer than 0.10 seconds



**Q-T Interval:** represents the time necessary for ventricular depolarization and repolarization.


- **Location:** Extends from the beginning of the QRS complex to the end of the T wave (includes the QRS complex, S-T segment, and the T wave)

**Duration:** Varies according to age, sex, and heart rate

**T Wave:** represents the depolarization of the ventricles. On rare occasions, a U wave can be seen following the T wave. The U wave reflects the depolarization of the His-Purkinje fibers.


- **Location:** Follows the S wave and the S-T segment  
**Amplitude:** 5mm or less in standard leads I, II, and III; 10mm or less in precordial leads V1-V6.

**Duration:** Not usually measured



**S-T Segment:** represents the end of the ventricular depolarization and the beginning of ventricular repolarization.

- **Location:** Extends from the end of the S wave to the beginning of the T wave  
**Duration:** Not usually measured



## QRS detector:

There are several applications which require an accurate detection of the QRS complex of the ECG. For example, arrhythmia monitors for ambulatory patients analyze the ECG in real time and when an arrhythmia occurs, the monitor stores a time segment of the abnormal ECG. This kind of monitor requires an accurate QRS recognition capability. Thus, QRS detection is an important part of many ECG signal processing systems.

Figure below show detection of QRS waves.





Fig1. QRS detection

Figure Below show the block Diagram of QRS Detector:



## Pulse missing detector:


When electrical excitation sometimes fails to pass through the A-V node or bundle of His, this intermittent occurrence is said to be called second degree heart block. Electrical conduction usually has a constant P-R interval, in the case of type 2 block atrial contractions are not regularly followed by ventricular contraction

Looking at the ECG you'll see that:

- Rhythm - Regular
- Rate - Normal or Slow
- QRS Duration - Prolonged
- P Wave - Ratio 2:1, 3:1
- P Wave rate - Normal but faster than QRS rate
- P-R Interval - Normal or prolonged but constant

Figure below show some missed QRS waves which can be detected by pulse missing detector which generate trigger pulse which indicate presence of arrhythmia and buzz the alarm.

:



## Conclusion:

## Results:

1. Observe the instance at which QRS is missing.

### ***Assignment:***

1. How do you detect QRS complex from ECG waveform ?
2. How do you detect various missing complexes from ECG waveform?

### ***Further Reading:***

1. Khandpur R. S., "Handbook of Biomedical Instrumentation", Tata McGraw-Hill, Second Edition.
2. Webster J. G., "Medical Instrumentation Application & Design", John Wiley & Sons, Inc., Fourth Edition.
3. Carr, J. J., Brown, J. M., 'Introduction to Biomedical Equipment Technology', Pearson Education, Fourth Edition.
4. Paul Horowitz, Winfield Hill., "The art of electronics", Cambridge University Press, 1989, Second Edition.
5. Ramakant A. Gayakwad, "Op-Amps and Linear Integrated Circuits," 4th edition, Pearson Education.